

## User's Checklist

**Serial Number:** \_\_\_\_\_

**AED Manufacturing Date:**

[illegible]

Status Indicator	
<b>Green:</b> the equipment is turned ON	Can work correctly
<b>Flashing green:</b> the equipment is in the standby status	Ready for operation at any time
<b>Flashing red and Beeping:</b> auto test failure is detected on the equipment.	<b>Defibrillator may still be operational</b> but still need to Immediately contact <a href="http://www.MindrayC1.com.au">www.MindrayC1.com.au</a>   1800 637 827   <a href="mailto:sales@MindrayC1.com.au">sales@MindrayC1.com.au</a>
<b>No Flashing Light and No Beeping:</b> no battery is installed or the battery is malfunctioning.	<b>DEFIBRILLATOR IS NOT OPERATIONAL.</b> Immediately contact <a href="http://www.MindrayC1.com.au">www.MindrayC1.com.au</a>   1800 637 827   <a href="mailto:sales@MindrayC1.com.au">sales@MindrayC1.com.au</a>

# Mindray® BeneHeart C Series Defibrillator

[illegible]

Initials	Print Full Name

Initials	Print Full Name